# JOVS

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# QUICK START GUIDE CARD

## RevivaEye Collagen Light Eye Patches

### [First use must do]

- Apply it on a small area behind the ear or inside the arm for 10 minutes first
- If there is no redness, stinging, or itching, then use it on the eye area.
- For the first use, strictly follow the steps and do not layer or extend the time.

### [Frequency of use]

- Use the 1-4 Day cycle set
- Once a week is sufficient. Do not frequently reuse the same area.

### [Safety Taboos] (Must-read)

- Pregnant women, breastfeeding women, during menstruation, and those under 18 years old are prohibited from using.
- Those wearing medical devices such as pacemakers / electrocardiographs / artificial heart-lung machines are prohibited from using.
- Those with eye inflammation, damage, post-operation, hypertension, heart disease, photosensitive skin are prohibited from using.
- Do not look directly at the light source; it is not waterproof, and it is strictly prohibited to immerse in water.
- Stay away from water sources, high temperatures, and corrosive environments; do not disassemble or repair by yourself.

### [Quick steps] \*This guideline provides instructions for the daily opening of the packaging. For more detailed usage methods, please refer to the instructions.

#### Use guide for DAY 1 /DAY 3

##### ① Hydration

**Use: Moisturizing Film**  
Apply the Moisturizing Film and leave it on for 5 minutes. Then remove it and dry the face.

##### ② Microcrystalline Effective Ingredients

**Use: Black Spot Essence Patch + Hyaluronic Acid Patch**  
First, fix the Black Spot Essence Patch on the under-eye area, then immediately cover the Hyaluronic Acid Patch. Gently press from the middle to both ends to make the attachment stable.  
\*Note: The Black Spot needs to penetrate the skin. Please do not touch the Black Spot with your fingers.

##### ③ Phototherapy

**Use: COLLAGEN LIGHT ENERGY CAPSULE**  
Before use, remove the transparent insulating paper under the battery and the white isolation paper under the liner, and stick it on the under-eye area. Keep it for 60-90 minutes to enjoy the treatment.

#### Use guide for DAY 2 / DAY 4

##### ① Keep moist

**Use: Essence Freeze-Dried Patch**  
Keep the area around the eyes slightly moist. Apply the Essence Freeze-Dried Patch to the area around the eyes.

##### ② Wet compress

**Use: Moisture-Enhancing Moisturizing Film**  
Apply the Moisture-Enhancing Moisturizing Film onto the Essence Freeze-Dried Patch and leave it on for 15 minutes. Then remove it.