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B l a c k e n P r o

【First use must do】

- First, test on a small area of the arm / behind the ear for 5-10 minutes.
- If there is no redness, pain, itching or other discomfort, then use it normally on the face / body.

【Frequency of use】

- Initial Adaptation: 2-3 times per week
- Regular Care: 3-5 times per week
- Avoid repeated use of the same area within one day.

【Safety Taboos】(Must-read)

- Pregnant women, people with photosensitive skin, and patients with skin diseases are prohibited from using.
- Those wearing medical electronic devices such as pacemakers are prohibited from using.
- Do not use on wounds, tattoos, or inflamed areas.
- Stay away from water sources, high temperatures, and corrosive environments.

【Mode】

BRI./PORE mode

- > Energy is gentle
- > Suitable for the entire face
- > Uniformly brightens skin tone and improves pores

REDNESS mode

- > Energy is moderate
- > Suitable for skin with corresponding characteristics (flushed skin)
- > Reduces the appearance of skin flushing

POWER mode

- > Energy is strong
- > Used for specific areas corresponding to the location of spots
- > Targeted at stubbornly hard-to-remove spots/pimples

【Quick steps】

- 1 Thoroughly clean the facial skin.
 - 2 Connect the power supply, press and hold the power button for 2 seconds to turn on the device, and briefly press the power button to set the desired function mode. Start using from the BRI./PORE mode.
 - 3 03 Please wear the special goggles correctly during operation.
 - 4 04 Align the light output port of the instrument with the position where light needs to be shone. When the sensing indicator light turns on, briefly press the light-shooting button to perform the light-shooting. (Press and hold the light-shooting button for 2 seconds to enter/exit the consecutive-shooting mode; if close to the skin around the eyes, it is recommended not to use consecutive-shooting, and perform single-shot light-shooting while wearing goggles and closing your eyes.)
 - 5 For areas with redness on the face, use the corresponding mode (for allergies or inflamed skin, it is prohibited); for areas with pigmentation, use the corresponding mode (for the position of moles, please use the mole-removing patch). ※ Choose the applicable mode based on your skin type.
- After the operation, apply a basic moisturizing and hydrating mask for 15-20 minutes.