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*Note: Avoid areas specified in the manual that cannot be exposed to light.

【Frequency of use】

Initial adaptation	2-3 times per week
Regular care	3-5 times per week

*Do not reuse the same area multiple times within one day.

【First use must do】

- Apply a small area behind your arm or ear for 5 to 10 minutes first
- If there is no redness, stinging, itching, etc., use it on the face / body.
- Gradually adapt throughout the process from the lowest setting.

【Safety Taboos】(Must-read)

- Pregnant women, breastfeeding women, and those with photosensitive skin are prohibited from using.
- Those wearing medical electronic devices such as pacemakers are prohibited from using.
- Do not use on wounds, tattoos, inflamed, or allergic areas.
- Stay away from water sources, high temperatures, humidity, and corrosive environments.

【Quick steps】


Thoroughly clean the facial skin

Connect the power supply and press the power-on button for 2 seconds to turn on. By default, it will enter the ACNE mode.


Place the light outlet of the instrument close to the area to be illuminated, then press the light-on button briefly to emit light once, double-click the light-on button to emit light continuously, and double-click again to exit the continuous flashing. *When illuminating, please close your eyes tightly to avoid looking directly at the light source.

Based on the condition of facial acne, press the power-on button briefly to switch to the corresponding mode. *If there is no operation for 3 consecutive minutes, the device will automatically shut down.

After the operation is completed, give the skin adequate hydration.

Keep physical sun protection for 48 hours after using this product. (It is not recommended to apply makeup when the skin is in an acne-prone state. If makeup is necessary, please use a light makeup.)